


Activated Charcoal filters are recommended for use with anaesthetic workstations:

1. During a suspected Malignant Hyperthermia crisis;
2. For patients undergoing anaesthesia who have a known or suspected increased risk of developing Malignant Hyperthermia.



**Guidelines based on evidence from:**

1. Bilmen JG, Gillies RL. *Anaesth Intensive Care* 2014; **42**: 51-58
2. Birgenheier N, Stoker R, Westonskow D, Orr J. *Anesth Analg* 2011; **112**: 1363-70.
3. Cottron N, Larcher C, Sommet A, Fesseau R, Alacoque X, Minville V, Fourcade O, Kern D. *Anesth Analg* 2014; **119**:67-75